

Field's Stain A & Field's Stain B with Fixative spray (Field's Stain Kit)**Intended use**

Ready to use combi kit each of Field's Stain A and Field's Stain B with Heme-Fixative for staining of blood smears.

Summary

Field's stain is a version of a Romanowsky stain, used for rapid processing of the specimens. Polychromated methylene blue and eosin stains selectively basophilic and acidophilic cellular elements to demonstrate blood cells and hemo-parasites. It is used for staining thick blood films to discover malarial parasites.

Principle

Field's Stains contain methylene blue and eosin. These basic and acidic dyes induce multiple colors when applied to cells. The fixative, methanol does not allow any further change in slide. The basic component of white cells (cytoplasm) is stained by acidic dye, and they are described as eosinophilic or acidophilic. The acidic component (nucleus with nuclei acid) takes blue to purple shades of the basic dye and is called basophilic. The neutral component of the cells is stained by both the dyes.

Reagents / Contents**Field's Stain A**

Methylene blue	1.45 g
Disodium hydrogen phosphate	9.09 g
Azur II	0.91 g
Potassium dihydrogen orthophosphate	11.36 g
Distilled water	977.2 mL

Field's Stain B

Eosin Y	2.60 g
Disodium hydrogen phosphate	10.0 g
Potassium dihydrogen orthophosphate	12.50 g
Distilled water	974.90 ml

**Formula adjusted, standardized to suit performance parameters

Appearance

Field's Stain A: Dark blue coloured solution

Field's Stain B: Orange coloured solution.

Storage and stability All reagents are stable at room temperature until the expiry date stated on each label. Do not spray fixative onto open fire. Avoid inhaling the spray.

Materials required but not provided

Clean grease-free glass slide, staining rack, blotting paper, immersion oil (Cat. No.207090110025), and microscope.

Type of Specimen: Clinical specimen: Blood samples.

Procedure

1. Prepare a thick blood film.
2. Fix blood film with Heme-Fixative for 1 minute and air dry.
3. Wash off fixative with tap water and air dry.
4. Cover the blood smear with Field's B stain and keep for 1 minute.
5. Gently wash the slide with tap water.
5. Immediately add an equal volume of Field's A stain and stain for 1 minute.
6. Rinse the slide with tap water and place upright to drain and dry.
7. Observe under microscope, 40X and 100X under oil immersion lens.

Interpretation of results

Neutrophil: Segmented polymorph with purple granules
Eosinophil: Polymorph with large orange granules
Basophil: Polymorph with large blue granules
Lymphocyte: Dark blue nuclei with clear blue cytoplasm
Monocyte: Kidney shaped nucleus with smoky blue cytoplasm
Platelet: Small pale bluish cells
Red blood cells: Pinkish red cells
Malaria Parasites: Pale bluish with Red Schuffner's Dots (Trophozoite)

Warranty

This product is designed to perform as described on the label and pack insert. The manufacturer disclaims any implied warranty of use and sale for any other purpose.

Reference

1. Data on file: Microxpress®, A Division of Tulip Diagnostics (P) Ltd.

Product Presentation:

Cat No.	Product	Pack Size
207060240250	Field's Stain Kit	2 x 250 mL
207060240500		2 x 500 mL
207060220250	Field's Stain A	250 mL
207060220500		500 mL
207060221000		1000 mL
207060225000		On request
207060230250	Field's Stain B	250 mL
207060230500		500 mL
207060231000		1000 mL
207060235000		5000 mL

Disclaimer

Information provided is based on our inhouse technical data on file, it is recommended that user should validate at his end for suitable use of the product.
