

Nyflock II™

Gamma-Irradiated Nylon microfibers flocked swab, recommended for nasopharyngeal and urethral collection with breakpoint.

Application

Nyflock II™ comprise an applicator shaft tip coated with nylon fibers that are arranged perpendicular fashion to optimize specimen collection and elution into transport media.

Technical Data Sheet

Parameters	UOM	Standard
Length	cm	15.5 ± 0.2
Distance from tip to break point	cm	10.0 ± 0.2
Weight of the swab	g	0.5 ± 0.02
Colour		Off White

MOC

Stick – Polypropylene
Swab - Nylon

Note – Gamma-Irradiation Certificate provided on request.

Reference

Data on file: Microexpress®, A Division of Tulip Diagnostics (P) Ltd.

Product Presentation:

Cat No.	Product	Pack Size
208140360100	Nyflock II™	100 Nos.

Disclaimer

Information provided is based on our inhouse technical data on file, it is recommended that user should validate at his end for suitable use of the product.
