

Poly-E swab I™

Gamma-Irradiated Polyester/Dacron® foam tipped swab, thermally bonded with polypropylene stick with breakpoint.

Application

Poly-E swab I™ Polyester/Dacron® swab is used for cleaning broad surfaces and flat areas and is widely recommended for surface sampling and cleaning validation work

Technical Data Sheet

Parameters	UOM	Standard
Head Length	cm	2.6 ± 1
Handle Length	cm	10.0 ± 1
Total Length	cm	12.6 ± 1
Weight of the swab	g	1.60 ± 0.05
Colour stick		Blue
Swab		White

MOC

Stick – Polypropylene

Swab - Dacron

Note – Gamma-Irradiation Certificate provided on request.

Reference

Data on file: Microxpress®, A Division of Tulip Diagnostics (P) Ltd.

Product Presentation:

Cat No.	Product	Pack Size
208160690100	Poly-E swab I™	100 Nos.

Disclaimer

Information provided is based on our inhouse technical data on file, it is recommended that user should validate at his end for suitable use of the product.
