

**Poly-E swab II™**

Gamma-Irradiated Polyester/Dacron® foam tipped swab, thermally bonded with polypropylene stick.

**Application**

Poly-E swab II™ Polyester/Dacron® swab is used for cleaning broad surfaces and flat areas and is widely recommended for surface sampling and cleaning validation work

**Technical Data Sheet**

Parameters	UOM	Standard
Head Length	cm	$1.7 \pm 0.1$
Handle Length	cm	$13.5 \pm 0.1$
Total Length	cm	$15.2 \pm 0.1$
Weight of the swab	g	$0.75 \pm 0.05$
Colour stick		Blue
Swab		White

**MOC**

Stick – Polypropylene

Swab - Dacron

**Note** – Gamma-Irradiation Certificate provided on request.

**Reference**

Data on file: Microxpress®, A Division of Tulip Diagnostics (P) Ltd.

**Product Presentation:**

Cat No.	Product	Pack Size
208160750100	Poly-E swab II™	100 Nos.

**Disclaimer**

Information provided is based on our inhouse technical data on file, it is recommended that user should validate at his end for suitable use of the product.

---