

**Steriloop - 10™**

Gamma-Irradiated sterile 10 µL inoculation loop and stabber tip.

**Application**

The loop is used in the cultivation of microbes on plates by streaking the inoculum on surface of the plate, whereas stabber tip is ideal for making stab inoculations into agar slants or tubes of solid culture medium.

**Technical Data Sheet**

Parameters	UOM	Standard
Length	cm	19.4 ± 0.2
Outer diameter of loop	cm	0.6 ± 0.05
Inner diameter of loop	cm	0.40 ± 0.02
Weight	g	0.95 ± 0.05
Colour		White

**MOC** - Polypropylene

**Note** – Gamma-Irradiation Certificate provided on request.

**Reference**

Data on file: Microxpress®, A Division of Tulip Diagnostics (P) Ltd.

**Product Presentation:**

Cat No.	Product	Pack Size
208191310100	Steriloop - 10™	100 Nos.

**Disclaimer**

Information provided is based on our inhouse technical data on file, it is recommended that user should validate at his end for suitable use of the product.

---