

**Sterispread™**

Gamma-Irradiated sterile disposable L- spreader.

**Application**

The L-shape of the spreader allows even application of samples/bacterial cultures on the entire petridish or plate surface.

**Technical Data Sheet**

Parameters	UOM	Standard
Length	cm	15.8 ± 0.2
Width	cm	5.7 ± 0.1
Thickness	cm	0.35 ± 0.02
Weight	g	2.6 ± 0.1
Colour		White

**MOC** - Polypropylene

**Note** – Gamma-Irradiation Certificate provided on request.

**Reference**

Data on file: Microxpress®, A Division of Tulip Diagnostics (P) Ltd.

**Product Presentation:**

Cat No.	Product	Pack Size
208191260100	Sterispread™	100 Nos.

**Disclaimer**

Information provided is based on our inhouse technical data on file, it is recommended that user should validate at his end for suitable use of the product.

---