

Steritrans-E™

Ethylene Oxide sterile polypropylene tube, cotton bud with polypropylene stick.

Application

Steritrans-E™ is used for collection and transport of specimens from the target area to Laboratory.

Technical Data Sheet**a) Transport Tube**

Parameters	UOM	Standard
Height	cm	14.8 ± 0.2
Outer diameter of mouth	cm	1.2 ± 0.1
Inner diameter of mouth	cm	1.1 ± 0.1
Outer body diameter of tube	cm	1.3 ± 0.1
Weight of the tube	g	3.0 ± 0.1
Colour		Transparent

MOC – Polypropylene tube

b) Yellow Colour Cap

Parameters	UOM	Standard
Height	cm	3.0 ± 0.1
Outer Body diameter	cm	1.3 ± 0.1
Weight of the cap	g	1.6 ± 0.1
Colour		Yellow

MOC – Polypropylene

c) Cotton Swab

Parameters	UOM	Standard
Length	cm	15.2 ± 0.1
Outer diameter of polypropylene stick	cm	0.24 ± 0.02
Weight of the swab	g	0.40 ± 0.1
Colour		White

MOC

Stick – Polypropylene
Swab - Cotton

Note - Gamma Irradiation Certificate provided on request.

Reference

Data on file: Microexpress®, A Division of Tulip Diagnostics (P) Ltd.

Product Presentation:

Cat No.	Product	Pack Size
208191290100	Steritrans-E™	100 Nos.

Disclaimer

Information provided is based on our inhouse technical data on file, it is recommended that user should validate at his end for suitable use of the product.
